


Attends Contours

(Shaped Pads)

Attends Contours are body contoured incontinence pads designed to sit close to the body. Advanced absorbent technology in the core ensures that all urine and unpleasant odours are locked quickly inside the pad. The soft flexible leakage barriers ensure maximum leakage protection. Each pad has a pale blue wetness indicator, which disappears when wet, to help you identify when the pad needs to be changed.

Attends Contours Fitting Guide




Apply the stretch pants, then pull them down to just above the knee and fold the waistband over.

Open the pad fully, making sure the back sheet is on the outside and the soft white lining is next to the skin.

Gently fold and cup the pad lengthways to ensure that the leakage barriers stand up. Don't touch the leakage barriers - they will cup and fit naturally with the pad.




Place the pad between the thighs if possible from front to back with the larger area of the pad towards the back of the body.



Smooth out the pad and pull up the stretch pants high into the groin to ensure a snug fit.



Adjust the pants to cover the pad and ensure that the back sheet is not folded inwards.



The pants should cover the pad. Finally check that both the pad and the stretch pants are fitted snugly.

Attends Contours dos and don'ts

- ✓ Assess the user for correct absorbency
- ✓ Measure the user around the waist for correct fixation pant size
- ✓ Always fold pad lengthways to create bowl shape
- ✓ Insert pad from front to back
- ✓ Ensure pad is fitting close to the body with wider part to the back
- ✓ Fixation pants should always be used with Comfort range
- ✗ Never place one pad inside another to increase absorbency
- ✗ Never use talc as it reduces absorbency of the product
- ✗ When creams are prescribed, use sparingly to ensure the product works effectively
- ✗ The products are designed to be body worn, do not use as bed or chair protection